

### All the cleaning/disinfecting power you need, right where you need it.

How did we improve the unparalleled cleaning/ disinfecting power and non-toxic safety of PathoSans? By making it more convenient than ever to use in daily work.

This new system is ideal for large facilities, offices and school campuses with multiple cleaning sites because you don't have to carry anything heavy. The 2.5 gallon transport containers are easily carried to any Dilution Station. Then, just make the amounts you need. No physical strain. No mess. No bulky containers to store.

- PS600-CT Generator produces PathoClean (cleaner/degreaser) and PathoCide (sanitizer/disinfectant) in concentrated formulas.
- Place PS600-CT Remote Dilution Stations at multiple sites where you do most of your cleaning.
- 3 Easily carry PathoClean and PathoCide concentrates to any Dilution Station in our lightweight 2.5 gallon containers and, within minutes, you have the world's best cleaning products ready-to-use.



# PATHOSANS • THE CLEAN THAT SAYS YOU CARE

### PS600-CT CONCENTRATE GENERATOR Remote Dilution Station

The PS600-CT Concentrate Generator produces concentrated formulas of PathoClean and PathoCide. You simply carry that concentrate in our lightweight 2.5 gallon transport containers to these Remote Dilution Stations and you're good-to-go. Place Dilution Stations everywhere you typically do a lot of cleaning. This chart shows how efficient this "cleaning ecosystem" will be in your large facility, offices or campus.

## Talk to us today about how PathoSans PS600-CT can work for you.

🕒 1.833.553.2648 🛛 🗖 PathoSans.com



#### PathoClean™


• One 2.5-gallon transport container of concentrate provides 25 gallons of ready-to-use PathoClean (9:1 dilution ratio).

#### PathoCide™

• One 2.5-gallon transport container of concentrate provides 17.5 gallons of ready-to-use PathoCide (6:1 dilution ratio).

#### **FEATURES**

- Generators come equipped with on-board 10 gallon tanks for both PathoClean and PathoCide.
- Each generator can accommodate up to four transport containers of each product at one time.
- One dilution station is included on each generator for easy access to ready-to-use chemistry directly from the unit.



# PATHOSANS • THE CLEAN THAT SAYS YOU CARE

